

# Patient Shield

Re-useable

## Spec Sheet



## Description

Specially designed to minimize the transfer of pathogens between patients and patient handlers, equipment and the environment. Requires the delivery of supplemental oxygen, preferably via non-rebreather, to prevent carbon dioxide build-up.

## Features

- Hook-and-loop fasteners secure the shield to the backrest frame to pull the shield around the sides of the cot
- Use with backrest raised or lowered
- Does not interfere with IV Poles or oxygen delivery

## Specifications

PART NUMBER	MANU 101255
COMPATIBILITY	Fits <b>ANY</b> North American Ambulance Cot
MATERIAL	Clear acrylic
CARE	Clean and disinfect according to local protocols. Replace if acrylic is torn or punctured.
<b>! IMPORTANT !</b>	Supplemental oxygen must be used to prevent carbon dioxide build-up.



Add: Via B. Zallone, 26  
Pieve di Cento - BO- 40066

Web: [www.ferno.it](http://www.ferno.it)

Ph: 051.6860028

Fax: 051.6861508