

Winter kit

cod.21-0120-019



 **FERNO**[®]


Via B. Zallone n. 26
40066 Pieve di Cento [Bo] | Italy
P +39 051 6860028 | F +39 051 6861508

info@ferno.it

www.ferno.it

Numero Verde
800 501711

 Winter Kit



 **FERNO**[®]

This is what we call
WINTER KIT.


The rigid shells you see here make the stretcher multifunctional and suitable also in the winter.

By adding the extension handles, the stretcher can be dragged on snow and therefore you can use it on the slopes and off-piste rescue operations.

The fin-like device under the stretcher makes it stable and manageable on snow.

The shells are made of fiberglass with a kevlar inner layer to increase resistance to abrasion.

- ▶ **EASY TO INSTALL WITH SIMPLE RUBBER BANDS.**
- ▶ **COMPATIBLE WITH DIVISIBLE TITAN-T STRETCHER.**
- ▶ **THREE PIECES DIVISIBLE HANDLES AND STOWED IN SMALL PLACE.**
- ▶ **SHELLS ARE EASY TO REPAIR IF BROKEN.**

 **Shell**



Fiberglass shell



▶ Winter kit

cod. 21-0120-019

no.1 head end shell, no.1 foot end shell,
handle extensions - 4pcs, no.1 connection bar

 **Handles**

