

40

Stair Chair

Code: 0714865





Lightweight, compact and sturdy design makes **model 40 stair chair** ideal for confined space and for a patient safe transport over most surfaces.

MAIN FEATURES

- Excellent transportability over most surfaces and excellent maneuverability in tight spaces thanks to two 127 mm stationary rear wheels.
- Equipped with two folding rear handles and two telescopic front handles. The front handles are adjustable to two positions for an ergonomic grip by operators.
- Vinyl-coated nylon cover of seat and backrest is simple to remove and resists stains, blood and bodily fluids guaranteeing an easy disinfection.
- Two quick-release patient restraints and leg restraint hold the patient securely.

SPECIFICATION

Height		910 mm
Width		510 mm
Depth	folded configuration	180 mm
	handles retracted	690 mm
	handles extended	1150 mm
Length (folded configuration)		910 mm
Weight		8,6 kg
Load Limit		159 kg

OPTIONAL ACCESSORIES

Description	Code
Fastener for model 40 stair chair	WS-40