

## Technical Data

Weight	1,6 kg	1,9 lbs
Width	300 mm	11,8 inch
Thickness	190 mm	7,5 inch
Height	165 mm	6,5 inch

## Lateral Blocks

Width	85 mm	3,3 inch
Thickness	230 mm	9,1 inch
Height	150 mm	5,9 inch

## Certification

- EN ISO 10993-17:2009 - *TESTED AND CERTIFIED material that does not release hazardous substances when in contact with skin and sweat.*
- EN ISO 10993-18:2020 - *Biocompatible material.*
- EN ISO 10993-5:2009 - *Non-cytotoxic material.*
- EN ISO 10993-23:2021 - *Non-irritating material for the skin.*
- COMPLIANT WITH THE NEW MDR REGULATION

## Accessories



Carrying and storage bag  
cod. 10-9900-005

 **FERNO®**

*Italy*

P +39 051 6860028  
[info.it@ferno.com](mailto:info.it@ferno.com)

[www.ferno.it](http://www.ferno.it)

*Schweiz*

P +41 41 259 60 00  
[info.ch@ferno.com](mailto:info.ch@ferno.com)

[www.ferno-schweiz.ch](http://www.ferno-schweiz.ch)

# B-lock2

 **FERNO®**

**-VELCRO**  
**+QUICK**  
**+CLEAN**  
**+SAFE**



 **FERNO®**

# B-lock2

FERNO



B-Soft, B-Safe, B-Lock!

## Head Immobilizer B-lock2

Lightweight, compact, and easy to apply. **Entirely manufactured in Italy**, the B-lock2 ensures complete and secure immobilization of the patient's head. When used in combination with the Ferno WizLoc cervical collar, it provides safe and comfortable transport for the patient.

The B-lock head immobilizer is equipped with two lateral blocks, a "K" strap and a "T" strap, allowing for quick immobilization.

The absence of velcro surfaces on both the base and lateral blocks provides optimal infection control and makes it easier to disinfect.

Enhanced comfort: The B-lock2 material is **SOFT TOUCH**, providing a better experience when in contact with the patients face as it adapts to its shape.

**TESTED:** The material is **100% recyclable** and does not release any hazardous substances when in contact with the skin.

The B-lock2 head immobilizer is made of a thermoplastic material with an **ANTIBACTERIAL ADDITIVE**, which prevents microbiological proliferation and improves hygienic storage conditions.

The head immobilizer base, in addition to providing stability to the lateral blocks, features an integrated V-shaped QHI that enhances comfort and ensures correct head positioning on the base.

The **ELLIPTICAL ERGONOMIC OPENINGS** on the lateral blocks of the B-lock make positioning the head immobilizer easier. This allows for direct contact with the patient, providing greater comfort around the ear area when subjected to immobilization pressure.



**Tested and biocompatible according to ISO 10993 standard.**

Using the new "K" strap and "T" strap, the B-lock2 can be secured to the stretcher for **PERFECT IMMOBILIZATION** of the patient.



Rev. 28012026

SIMPLE AND OPTIMAL SANITIZATION

FERNO