

# FLOTATION COLLAR

*Water Rescue*

Code: 11-0124A



# FLOTATION COLLAR

**FERNO**  
When It's Critical®



The lightweight Flotation Collar provides flotation assistance during patient handling in water.

It is flexible, two piece design allows for quick and easy installation on a wide variety of basket stretchers.

## MAIN FEATURES

- *Designed for use with basket stretchers.*
- *45 lb. of bouyancy to keep patient afloat.*
- *Conveniently stores inside basket stretcher.*
- *Compatible with Model 71, Gazelle Series, Titan Series, Advantage and Spartan.*

## SPECIFICATIONS (PER COLLAR)

Length	1300 mm
Width	100 mm
Depth	200 mm
Weight	1,4 kg

