

Industry Standard Extrication





KENDRICK EXTRICATION DEVICE

The Original Ferno K.E.D.® is a versatile, improved means of immobilizing and extricating patients from auto accidents or confined spaces. The wrap-around design provides horizontal flexibility for easy application and vertical rigidity for maximum support of the spine, neck, and head during extrication. The K.E.D.® immobilizes the torso, head and neck, enabling prompt extrication while minimizing risks of further injury. X-rays or advanced life-support procedures can be performed with the K.E.D.® in place.

The Ferno K.E.D.® includes wrap-around vest, Adjusta-Pad neck roll, two head straps and a carrying case. Heavy-duty, vinyl-coated nylon is durable and easy to clean and the K.E.D.® stores rolled up, ready for immediate use.



FEATURES

- Proven means of immobilizing and extricating patients from auto accidents or confined spaces
- The wrap-around design provides vertical rigidity for maximum support when immobilizing
- the spine, neck, and head during extrication, minimizing the risk of further injury
- X-ray translucent
- Heavy-duty, vinyl-coated nylon is durable and easy to clean
- Case is enlarged to accommodate extrication collars
- Built-in handles enable rescuers to get a firm grip on the patient and K.E.D.[®]
 as one unit
- Adapts to a hip and pelvic splint by simple inversion
- Adjustable fold-back sides permit easy access to patient's chest
- Can be used for multiple patient sizes, including children and pregnant women

SPECIFICATIONS		
Length:	Open:	830 mm
	Rolled:	900 mm
Width:	Open:	800 mm
	Rolled:	270 mm
Thickness:	Open:	30 mm
	Rolled:	140 mm
Weight:		3 kg
Load Limit:		227 kg



